

SUGCT Antibody



PACO08114

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SUGCT

Uniprot

Q9HAC7

Synonyms:

chromosome 7 open reading frame 10; C7orf10; DERP13; FLJ11808; FLJ53872; ORF19 ;

Immunogen:

Human C7orf10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A