

# CPED1 Antibody



PACO08116

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CPED1

**Uniprot**

A4D0V7

**Synonyms:**

chromosome 7 open reading frame 58; C7orf58; FLJ21986; FLJ26813 ;

**Immunogen:**

Human C7orf58.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A