

C9orf78 Antibody



PACO08118

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

C9orf78

Uniprot

Q9NZ63

Synonyms:

chromosome 9 open reading frame 78; C9orf78; HSPC220; bA409K20.3 ;

Immunogen:

Human C9orf78.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A