

CACNB1 Antibody



PACO08135

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CACNB1

Uniprot

Q02641

Synonyms:

calcium channel, voltage-dependent, beta 1 subunit; CACNB1; CAB1; CACNLB1; CCHLB1; MGC41896 ;

Immunogen:

Human CACNB1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A