

CACNG5 Antibody



PACO08137

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CACNG5

Uniprot

Q9UF02

Synonyms:

calcium channel, voltage-dependent, gamma subunit 5; CACNG5; MGC126656; MGC126682 ;

Immunogen:

Human CACNG5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A