

CACYBP Antibody



PACO08138

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CACYBP

Uniprot

Q9HB71

Synonyms:

calcyclin binding protein; CACYBP; GIG5; MGC87971; PNAS-107; RP1-102G20.6; S100A6BP; SIP ;

Immunogen:

Human CACYBP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A