

CALR3 Antibody



PACO08151

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CALR3

Uniprot

Q96L12

Synonyms:

calreticulin 3; CALR3; CRT2; FLJ25355; MGC26577

Immunogen:

Human CALR3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A