

CAMP Antibody



PACO08161

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CAMP

Uniprot

P49913

Synonyms:

cathelicidin antimicrobial peptide; CAMP; CAP18; CRAMP; FALL-39; FALL39; HSD26; LL37 ;

Immunogen:

Human CAMP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A