

CAPRIN2 Antibody



PACO08178

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

CAPRIN2

Uniprot

Q6IMN6

Synonyms:

caprin family member 2; CAPRIN2; C1QDC1; EEG-1; EEG1; FLJ11391; FLJ22569; KIAA1873; MGC102894; MGC134847; MGC134848; caprin-2 ;

Immunogen:

Human CAPRIN2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A