

# CASP4 Antibody



PACO08197

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CASP4

**Uniprot**

P49662

**Synonyms:**

caspase 4, apoptosis-related cysteine peptidase; CASP4; ICE(rel)II; ICEREL-II; ICH-2; Mih1/TX; TX ; caspase4

**Immunogen:**

Human CASP4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A