

CCDC42 Antibody



PACO08228

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CCDC42

Uniprot

Q96M95

Synonyms:

coiled-coil domain containing 42; CCDC42; FLJ32734 ;

Immunogen:

Human CCDC42.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A