

CCDC51 Antibody



PACO08231

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CCDC51

Uniprot

Q96ER9

Synonyms:

coiled-coil domain containing 51; CCDC51; FLJ12436 ;

Immunogen:

Human CCDC51.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A