

PACO08234

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CCDC86

Uniprot

Q9H6F5

Synonyms:

coiled-coil domain containing 86; CCDC86; FLJ22321; MGC2574 ;

Immunogen:

Human CCDC86.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A