

CCRL2 Antibody



PACO08260

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CCRL2

Uniprot

O00421

Synonyms:

chemokine (C-C motif) receptor-like 2; CCRL2; CKRX; CRAM-A; CRAM-B; FLJ55815; HCR; MGC116710; MGC34104 ;

Immunogen:

Human CCRL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A