

CCT6B Antibody



PACO08264

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CCT6B

Uniprot

Q92526

Synonyms:

chaperonin containing TCP1, subunit 6B (zeta 2); CCT6B; CCT-zeta-2; CCTZ-2; Cctz2; TCP-1-zeta-2; TSA303 ;

Immunogen:

Human CCT6B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A