

# CD151 Antibody



PACO08269

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CD151

**Uniprot**

P48509

**Synonyms:**

CD151 molecule (Raph blood group); CD151; GP27; MER2; PETA-3; RAPH; SFA1; TSPAN24 ;

**Immunogen:**

Human CD151.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A