

CD164 Antibody



PACO08270

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CD164

Uniprot

Q04900

Synonyms:

CD164 molecule, sialomucin; CD164; MGC-24; MUC-24; endolyn ;

Immunogen:

Human CD164.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A