

CD24 Antibody



PACO08275

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CD24

Uniprot

P25063

Synonyms:

CD24 molecule; CD24; CD24A; FLJ22950; FLJ43543; MGC75043 ;

Immunogen:

Human CD24.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A