

# CD276 Antibody



PACO08280

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

CD276

**Uniprot**

Q5ZPR3

**Synonyms:**

CD276 molecule; CD276; B7-H3; B7H3 ;

**Immunogen:**

Human CD276.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A