

CD300LG Antibody



PACO08282

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CD300LG

Uniprot

Q6UXG3

Synonyms:

CD300 molecule-like family member g; CD300LG; CLM9; NEPMUCIN; TREM4 ;

Immunogen:

Human CD300LG.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A