

CD40 Antibody



PACO08290

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CD40

Uniprot

P25942

Synonyms:

CD40 molecule, TNF receptor superfamily member 5; CD40; Bp50; CDW40; MGC9013; TNFRSF5; p50 ;

Immunogen:

Human CD40.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A