

CD47 Antibody



PACO08294

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CD47

Uniprot

Q08722

Synonyms:

CD47 molecule; CD47; IAP; MER6; OA3 ;

Immunogen:

Human CD47.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A