

CD58 Antibody



PACO08297

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CD58

Uniprot

P19256

Synonyms:

CD58 molecule; CD58; FLJ23181; FLJ43722; LFA-3; LFA3 ;

Immunogen:

Human CD58.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A