

CD86 Antibody



PACO08303

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CD86

Uniprot

P42081

Synonyms:

CD86 molecule; CD86; B7-2; B70; CD28LG2; LAB72; MGC34413 ;

Immunogen:

Human CD86.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A