

CDK19 Antibody



PACO08314

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CDK19

Uniprot

Q9BWU1

Synonyms:

cell division cycle 2-like 6 (CDK8-like); CDC2L6; CDK11; KIAA1028; bA346C16.3 ; CDK19;

Immunogen:

Human CDC2L6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A