

# CDC34 Antibody



PACO08315

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDC34

**Uniprot**

P49427

**Synonyms:**

cell division cycle 34 homolog (S. cerevisiae); CDC34; E2-CDC34; UBC3; UBE2R1 ;

**Immunogen:**

Human CDC34.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A