

CDC37 Antibody



PACO08316

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CDC37

Uniprot

Q16543

Synonyms:

cell division cycle 37 homolog (S. cerevisiae); CDC37; P50CDC37 ;

Immunogen:

Human CDC37.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A