

CDC42 Antibody



PACO08318

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CDC42

Uniprot

P60953

Synonyms:

cell division cycle 42 (GTP binding protein, 25kDa); CDC42; CDC42Hs; G25K ;

Immunogen:

Human CDC42.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A