

CDC42EP4 Antibody



PACO08321

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CDC42EP4

Uniprot

Q9H3Q1

Synonyms:

CDC42 effector protein (Rho GTPase binding) 4; CDC42EP4; BORG4; CEP4; KAIA1777; MGC17125; MGC3740 ;

Immunogen:

Human CDC42EP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A