

CDC45 Antibody



PACO08322

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CDC45

Uniprot

O75419

Synonyms:

CDC45 cell division cycle 45-like (S. cerevisiae); CDC45L; CDC45; CDC45L2; PORC-PI-1 ;

Immunogen:

Human CDC45L.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A