

PACO08324

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

CDC6

Uniprot

Q99741

Synonyms:

cell division cycle 6 homolog (S. cerevisiae); CDC6; CDC18L; HsCDC18; HsCDC6 ;

Immunogen:

Human CDC6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A