

CDH13 Antibody



PACO08329

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CDH13

Uniprot

P55290

Synonyms:

cadherin 13, H-cadherin (heart); CDH13; CDHH ; T-cadherin, T-cad, P105, H-cadherin

Immunogen:

Human CDH13.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A