

# CDK2AP1 Antibody



PACO08336

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

CDK2AP1

**Uniprot**

O14519

**Synonyms:**

cyclin-dependent kinase 2 associated protein 1; CDK2AP1; DOC1; DORC1; ST19; doc-1; p12DOC-1 ;

**Immunogen:**

Human CDK2AP1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A