

CDK3 Antibody



PACO08337

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CDK3

Uniprot

Q00526

Synonyms:

cyclin-dependent kinase 3;

Immunogen:

Human CDK3.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A