

# CDK5RAP3 Antibody



PACO08341

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CDK5RAP3

**Uniprot**

Q96JB5

**Synonyms:**

CDK5 regulatory subunit associated protein 3; CDK5RAP3; C53; HSF-27; IC53; LZAP; MST016; OK/SW-cl.114 ;

**Immunogen:**

Human CDK5RAP3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A