

CDK9 Antibody



PACO08344

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CDK9

Uniprot

P50750

Synonyms:

cyclin-dependent kinase 9; CDK9; C-2k; CDC2L4; CTK1; PITALRE; TAK ;

Immunogen:

Human CDK9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A