

CDKN2B Antibody



PACO08348

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CDKN2B

Uniprot

P42772

Synonyms:

cyclin-dependent kinase inhibitor 2B (p15, inhibits CDK4); CDKN2B; CDK4I; INK4B; MTS2; P15; TP15; p15INK4b ;

Immunogen:

Human CDKN2B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A