

CEACAM21 Antibody



PACO08358

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

CEACAM21

Uniprot

Q3KPI0

Synonyms:

carcinoembryonic antigen-related cell adhesion molecule 21; CEACAM21; CEACAM3; FLJ13540; MGC119874; R29124_1 ;

Immunogen:

Human CEACAM21.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A