

CENPJ Antibody



PACO08364

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:**Protein Background:****Gene ID:**

CENPJ

Uniprot

Q9HC77

Synonyms:

centromere protein J; CENPJ; BM032; CPAP; LAP; LIP1; MCPH6; MGC131581; MGC131582; MGC142222; MGC142224 ;

Immunogen:

Human CENPJ.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A