

CEP63 Antibody



PACO08367

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

CEP63

Uniprot

Q96MT8

Synonyms:

centrosomal protein 63kDa; CEP63; FLJ13386; MGC78416 ;

Immunogen:

Human CEP63.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A