

# CEP68 Antibody



PACO08368

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

CEP68

**Uniprot**

Q76N32

**Synonyms:**

centrosomal protein 68kDa; CEP68; FLJ25920; FLJ36750; KIAA0582 ;

**Immunogen:**

Human CEP68.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A