
Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CFB

Uniprot

P00751

Synonyms:

complement factor B; CFB; BF; BFD; CFAB; FB; FBI12; GBG; H2-Bf; PBF2 ;

Immunogen:

Human CFB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A