

CHCHD7 Antibody



PACO08396

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CHCHD7

Uniprot

Q9BUK0

Synonyms:

coiled-coil-helix-coiled-coil-helix domain containing 7; CHCHD7; FLJ40966; MGC2217 ;

Immunogen:

Human CHCHD7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A