

CHEK2 Antibody



PACO08402

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CHEK2

Uniprot

O96017

Synonyms:

CHK2 checkpoint homolog (S. pombe); CHEK2; CDS1; CHK2; HuCds1; LFS2; PP1425; RAD53 ;

Immunogen:

Human CHEK2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A