

CIB4 Antibody



PACO08437

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CIB4

Uniprot

A0PJX0

Synonyms:

calcium and integrin binding family member 4; CIB4; KIP4 ;

Immunogen:

Human CIB4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A