

CKS2 Antibody



PACO08449

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

CKS2

Uniprot

P33552

Synonyms:

CDC28 protein kinase regulatory subunit 2; CKS2; CKSHS2 ;

Immunogen:

Human CKS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A