

CLCN1 Antibody



PACO08453

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

CLCN1

Uniprot

P35523

Synonyms:

chloride channel 1, skeletal muscle; CLCN1; CLC1; MGC138361; MGC142055 ;

Immunogen:

Human CLCN1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A