

CLCNKA Antibody



PACO08454

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CLCNKA

Uniprot

P51800

Synonyms:

chloride channel Ka; CLCNKA; CLCK1; CIC-K1; MGC61490; hCIC-Ka ;

Immunogen:

Human CLCNKA.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A