

# CLEC18A Antibody



PACO08460

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CLEC18A

**Uniprot**

A5D8T8

**Synonyms:**

C-type lectin domain family 18, member A; CLEC18A; FLJ77955; MGC164804; MRCL; MRCL1 ; MRLP2

**Immunogen:**

Human CLEC18A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A