

CLEC1A Antibody



PACO08461

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

CLEC1A

Uniprot

Q8NC01

Synonyms:

C-type lectin domain family 1, member A; CLEC1A; CLEC1; MGC34328 ;

Immunogen:

Human CLEC1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A